
# VIVASIGHT

By ETVIEW

Instructions for Use

## **Left-sided VivaSight-DL**

For use by trained personnel only

---

**For use only with Ambu® aView™ monitor**

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## VIVASIGHT-DL

**STERILE (ETO) unless package has been opened or damaged.**

### DO NOT REUSE.

Federal (US) law restricts this device to sale by or on the order of a physician.

### DESCRIPTION

VivaSight-DL is a single-use, left-sided, double-lumen cuffed PVC endobronchial tube with an embedded video imaging device and light source at the distal end of the tracheal lumen, and an integrated single use video/power cable (adapter cable) with connector. VivaSight-DL displays images of the airway onto the Ambu® aView™ monitor for as long as the device remains in place during intubation.

VivaSight-DL has two color-coded low-pressure cuffs with corresponding pilot balloons: clear cuff and balloon for tracheal blocking; blue cuff and balloon for bronchial blocking. Additionally, for cleaning the imaging lens, the VivaSight-DL has an injection port leading to two lumens along the tube's wall, which opens at the distal end of the imaging lens.

### INDICATIONS FOR USE

The ETVivaSight-DL System is intended to isolate the left or right lung of a patient for intensive care or surgery, one lung ventilation or one lung anesthesia.

The VivaSight-DL system is indicated for verifying tube placement and repositioning.

### CONTRAINDICATIONS

Do not use laser equipment in the immediate proximity of the VivaSight-DL. Do not use VivaSight-DL for patients with obstructions or stenosis in the left main bronchus.

## ACCESSORIES

Styler.

### WARNINGS

- VivaSight-DL is a single use product. Do not reuse: After use, handle and dispose according to local regulations for disposal of infected medical devices. Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.
- Do not attempt to resterilize any part of VivaSight-DL.
- Do not connect VivaSight-DL into any device or electrical outlet other than those supplied by ETViva or Ambu.
- The intended performance of the VivaSight-DL is reached only with the disposable devices and accessories provided or specified by Ambu or ETViva.
- VivaSight-DL is available in tube sizes of 35 Fr; 37 Fr; 39 Fr and 41 Fr (outer diameter - OD) and should be used for patients who require one of these sizes.
- Do not use the injection port for suctioning the patient.
- The cuff pressure should normally not exceed 25 cm H<sub>2</sub>O. Diffusion of a nitrous oxide mixture, oxygen, or air may either increase or decrease the cuff pressure.
- Do not use the camera in the event that inspired gas temperature exceeds 34 °C for an artificially ventilated patient.
- Do not overinflate the cuff. Overinflation can result in tracheal/bronchial damage, rupture of the cuff with subsequent deflation, or cuff distortion, which may lead to airway blockage.
- Do not use the tube if the cuff is damaged. Care must be taken

to avoid damaging the cuff during intubation.

- Do not use the tube if the cuff is damaged. Various bony anatomical structures such as teeth, or any intubation aid with sharp surfaces in the intubation path can damage cuff integrity. Care must be taken to avoid damaging the thin-walled cuff during intubation, which may result in the patient requiring traumatic extubation and reintubation.
- Both cuffs need to be completely deflated prior to repositioning the tube. Movement of the tube with inflated cuffs may result in damage to the cuff or in patient injury, which may result in medical interventions.
- If the tube is lubricated prior to intubation, it is essential to verify that lubricant has not entered or occluded the tube lumen or cuff-filling system, thereby preventing ventilation or damaging the cuff. Ventilation and proper functioning of the cuff may be impaired.
- Reliance on the graduated black depth marks on the tube should not substitute for expert judgment. The user should be aware of anatomical variations, including the length of the airway. Intubation and extubation should be performed following currently accepted medical techniques.
- VivaSight-DL's location should be verified every time the patient is moved. Should extreme flexion of the head (chin-to-chest) or movement of the patient (e.g., to lateral or prone positions) occur after intubation, ensure that VivaSight-DL remains in place.
- VivaSight-DL is for use in non-MRI environment only.
- Do not use an intubation stylet other than the stylet provided with the VivaSight-DL.
- To ensure ventilation as intended, be careful to fully align the arrows on the Y-connector lumen and rotator pointing either in the same direction or opposite direction depending on intended ventilation.
- Electronic equipment and the VivaSight-DL system may affect the normal function of each other. If the VivaSight-DL system is used adjacent to or stacked with other equipment, observe and verify normal operation of both the VivaSight-DL system and the other electronic equipment prior to using it. It may be necessary to adopt procedures for mitigation, such as re-orientation or relocation of the equipment or shielding of the room in which it is used. Consult the tables in appendix 1 (English version) for guidance in placing the VivaSight-DL system.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- The VivaSight-DL is not to be used when delivering flammable anaesthetics to the patient. This could potentially cause patient injury.
- Be careful to check whether the image on the screen is a live image or a recorded image.
- The surface temperature of the VivaSight-DL is likely to reach above 43 °C when the tube is placed outside the patient. Therefore switch off the monitor after the preuse test and switch it back on immediately before use, to prevent the risk of having an impact on the tissue. The surface temperature on the Adapter Cable is likely to reach between 41 °C and 43 °C.
- Prior to use, check that the pouch seal is intact and there are no impurities or damage on the product such as rough surfaces, sharp edges or protrusions which may harm the patient.

**CAUTIONS**

- Use of HF equipment in the immediate proximity of the VivaSight-DL may lead to impaired image quality.
- Adapter Cable and connectors shall be kept away from patient face to avoid facial harm.
- VivaSight-DL is intended for use by trained personnel only.
- Use the same size VivaSight-DL as the required standard tube size. Expert clinical judgment should be used in choosing the suitable tracheal tube size for each patient.
- Insert the stylet from the bronchial side of the tube only.
- Before use, always check for compatibility between the VivaSight-DL and the external instruments (e.g. bronchoscope, and suction catheter).
- Do not immerse the tube in liquid.
- The use of lidocaine topical aerosols has been associated with the formation of pinholes in PVC cuffs\*. To prevent cuff leaks, expert clinical judgment is required when using lidocaine topical aerosols.
- Use only water-soluble lubricants. Other lubricants might affect the cuffs.
- Avoid applying lubricants to the front of the camera or the inside of the tube, due to risk of affecting ventilation.
- Use a cuff pressure gauge to help monitor and adjust the cuff pressure. Inflation of the cuff by "feel" alone or by using a measured volume of air is not recommended as resistance is an unreliable guide during inflation.
- Syringes, stopcocks or other devices should not be left in the inflation system for an extended period.
- When connecting the Y-Connector, there is no need to remove the standard connectors from VivaSight-DL.

- To ensure secure placement once the connector or the Y-connector has been loosened or removed, the connector or the Y-connector and the corresponding part of the tube shall be cleaned with ethanol prior to reconnecting the connector or the Y-connector.
- Use of a lubricant solution to ease the insertion of the 15 mm connectors is not recommended as it may contribute to a accidental disconnection.
- Do not connect VivaSight-DL to any external monitor equipment other than aView.

\*Jayasuriya KD, Watson WF: "P.V.C. cuffs and Lidocaine-based aerosol"; Br J Anaesth. 1981 Dec; 53 (12): 1368.


**COMPATIBILITY**

When VivaSight-DL is used with external instruments (e.g. bronchoscope, suction catheter), use the following table of VivaSight-DL effective inside diameter (defined by ISO 16628) measurements as guidance for choosing the appropriate diameter of the external instrument that will be inserted through VivaSight-DL.

PRODUCT	Effective Inside Diameter
VivaSight-DL 35 Fr	4.0 mm
VivaSight-DL 37 Fr	4.2 mm
VivaSight-DL 39 Fr	4.8 mm
VivaSight-DL 41 Fr	5.0 mm

**Y-CONNECTOR OPERATION**

Make sure to connect the blue tube of the Y-connector to the blue bronchial tube of the VivaSight-DL and the clear tube of the Y-connector to the clear tracheal tube of the VivaSight-DL.

- The Y-connector has an arrow printed on both the tracheal and the bronchial lumen indicating the flow of air.
- An identical arrow is printed on the rotator connected to each of the two lumens.
- When the rotator arrow is pointing in the same direction as the arrow printed on the lumen, the lumen is open for ventilation.
- In order to close the lumen for ventilation, turn the rotator 180° degrees until the arrows are pointing in the opposite direction.

#### STORAGE AND TRANSPORTATION

- Store and transport VivaSight-DL at temperatures between 0 °C and 42 °C, relative humidity between 10-100% and atmospheric pressure between 80-109 kPa.
- Store in a dry, cool and dark place.

#### ADVERSE EVENTS

Adverse events associated with the use of VivaSight-DL are the same as those of standard double-lumen endobronchial tubes. The most common are: laryngospasm; vocal cord paralysis; injuries to lips, gums, tongue, teeth; and aspiration of gastric contents. Consult scientific literature for specific adverse reaction information.

Adverse events associated with use of standard double-lumen endobronchial tubes, are fracture or dislocation of the cervical spine, endobronchial or esophageal intubation, perforation of the trachea or esophagus, failed intubation, tube dislocation and misplacement of the bronchial tube.

#### TECHNICAL SPECIFICATIONS

The VivaSight-DL is an electrical safety defibrillation-proof type BF applied part.

- Imaging sensor: CMOS
- Resolution: CIF 320 x 240
- Video format: NTSC Composite Video Baseband Signal (CVBS)
- Light source: 2 LEDs (integrated)
- FOV: Horizontal ~76°, Vertical ~56°, Diagonal ~100°
- Focusing range: 12 mm-60 mm
- Power: up to 36 mA @ DC 5V
- Push VivaSight-EB into VivaSight-SL or equivalent.
- Operating ambient temperature: 10-37 °C (50-98 °F)
- Operating relative humidity: 30-75%
- Operating atmospheric pressure: 80-109 kPa
- Operating altitude: ≤2000 m

#### INTUBATION RECOMMENDATIONS

Intubation with VivaSight-DL is performed according to currently accepted medical techniques.

The tube position should be inside the trachea and inserted into the left main bronchus.

- Laryngoscope: Use of laryngoscope for intubation of VivaSight-DL.
- Stylet: The accompanying stylet is used to make VivaSight-DL more rigid, for use with or without laryngoscope.
- Oral: Approved for oral intubation procedures.

#### PRE-INTUBATION PREPARATION CHOOSING VIVASIGHT-DL SIZE

- Use the same size VivaSight-DL as standard required double-lumen endobronchial tube size.
- Use expert clinical judgment when choosing tracheal tube size for each patient.


## CONNECTING THE TUBE, VERIFICATION STEPS & PREPARATION FOR INTUBATION

1. Connect VivaSight-DL to the single use adapter cable. Be careful to align the white arrow on the single use adapter cable with the arrow on the VivaSight-DL tube round connector.
2. Connect the single use adapter cable to the aView™ monitor, be careful to align the arrow on the cable with the arrow indicator on the aView™ monitor. Verify that all electrical connections are stable and secure. Refer to aView™ monitor instructions for use.
3. Press the power button on top of the aView™ monitor for at least one second until a live image appears. Refer to Ambu® aView™ monitor instructions for use. The aView™ user interface will be ready after approximately one minute, when the hourglass symbol disappears.
4. Check that the two LEDs at the tip of the tube are illuminated. If they are not illuminated, replace the tube with another.
5. Verify that there is a stable image on the aView™ monitor (projected from the tip of VivaSight-DL). Then turn off the monitor again until immediately prior to intubation.
6. Test the cuffs for integrity by inflating and deflating them completely.
7. Become familiar with the feel of the blue and clear pilot balloons.
8. If applicable, verify that external instruments of appropriate size can be passed through the VivaSight-DL without resistance. There is no guarantee that instruments selected solely using the effective inside diameter will be compatible in combination with the VivaSight-DL.
9. Induce anesthesia, as per protocol.
10. Suction the oral cavity and airway area before inserting VivaSight-DL.


11. Prepare the patient, and lubricate the tube (if needed), avoiding the area around the camera lens at the distal end of VivaSight-DL.
12. Intubate the patient and remove the stylet. In case of obstructed or obscured view, see “Cleaning the Lens”.

## INTUBATION TECHNIQUE WITH VIVASIGHT-DL

1. Advance the VivaSight-DL bronchial tube through the vocal cords with the tip facing upward (figure 1).


2. At the proximal portion of the tube, the bronchial arm of the tube will be positioned above the tracheal arm (figure 2).


3. When the tip passes the vocal cords, and before the camera passes it, turn the tube 90 degrees counterclockwise until the camera on the tube is facing up pointing towards the vocal cords (figure 3). In this position, the tube passes through the vocal cords with the camera on the superior portion of the tube (figure 3).

4. When the camera passes the vocal cords, advance the tube further until final positioning, without additional tube maneuvering. When the tube is in place, the proximal portion of the bronchial arm will be horizontal and on the left side (figure 4).

When required, use a video laryngoscope to assist with correct intubation with VivaSight-DL.

## DEPTH SCALING

*Note: Suggested depth is 2 cm to 3 cm above the main carina. Keep in mind that VivaSight-DL's electronic wire is outside the protective lumen at 28 cm to 30 cm.*

Like with most intubation tubes, a series of depth marks appear on the side of VivaSight-DL.

Depth marks range from 25 cm to 31 cm, in 2 cm increments.

To determine a specific depth, gauge from the carina:

1. Insert VivaSight-DL's camera up to the depth of the carina.
2. Look at VivaSight-DL's depth marks on the side of the tube.

3. Pull VivaSight-DL back from the carina to the desired depth.

## VISIBILITY DURING PROCEDURE

In the unlikely event that the image is not shown during a procedure, disconnect the VivaSight-DL from the single use adapter cable to the aView™ monitor and continue to use VivaSight-DL as a standard airway tube. **Do not attempt to fix the connection or replace the monitor.** Turn off the aView™ monitor and dispose of the single use adapter cable.

## CLEANING THE LENS

If VivaSight-DL's lens becomes soiled or obscured by secretions, it can be cleaned by injecting air, prescribed liquid medications, or saline, where permitted by institutional policy, into the injection port of the rinsing tube.

Recommended lens cleaning method:

1. Inject 20 cc air into the injection tube port (RED), and then check image clarity.
2. If the image is not sufficiently clear, connect 10 cc syringe filled with 2 cc of Saline, to the flushing port and push the plunger. Perform this step twice (total of 4 cc of Saline).
3. Fill a 10 cc syringe with 5-10 cc of air, connect it to the flushing port and push the plunger. Perform this step twice.
4. Connect an empty 10 cc syringe, to the flushing port and suck the Saline. Perform this step twice.

## ADDITIONAL NOTES

- Continuous viewing of the airway enables real-time verification of the tube and the bronchial cuff position, early detection of adverse airway events, and secretion management under visual guidance.
- Validation: During static situations where VivaSight-DL is in place, and the patient is not being moved, validation that the system

is operating must be performed periodically by observing breathing movements and cardiac pumping movement in the image on the aView™ monitor.

- The device usage time function on the aView™ monitor is not available for VivaSight-DL. For this device, the device usage time will always display "00:00:00". This will not affect the functioning of the system.

## EXTUBATION

- Deflate the tracheal and endobronchial cuffs completely, as per protocol.
- Extubation is performed slowly, while examining the airway on the aView™ monitor one last time (document or record, as necessary).

## DISCONNECTING VIVASIGHT-DL

1. Disconnect the single use adapter cable from the aView™ monitor.
2. Turn off the aView™ monitor by pressing the power button for at least two seconds.

## DISPOSAL

Discard the single use VivaSight-DL tube, Y-connector, stylet and adapter cable in a suitable biohazard receptacle in a accordance with local regulations.

For disposal instructions for the aView™ monitor, refer to aView™ monitor instructions for use.

## MANUFACTURED BY:

ETView Ltd.

Catom 2 Street

Misgav Business Park

M.P. Misgav 2017900, Israel

www.etview.com

LB10097 REV03

 0483	
 Manufacture Date	 Expiry Date
	Consult instructions for use
	Do not use if package is damaged
	Do not reuse
<b>STERILE</b> <b>EO</b>	Sterilized Using Ethylene Oxide
	Non-sterile (single use adapter cable only)
<b>Rx ONLY</b>	US Federal law restricts this device to sale by or on the order of a physician
	Defibrillation-proof type BF applied part
	This product is not made with phthalates
	Temperature for transportation and storage 0°C / 32°F to 42°C / 107°F
	Humidity for transportation and storage 10% to 100%
	Pressure for transportation and storage 80kPa to 100kPa
	 Manufacturer
<b>EC REP</b>	Authorized Representative in the European Community
	Waste Bin Symbol
	This product is not made with natural rubber latex
<b>REF</b>	Reference Number
<b>LOT</b>	Lot Number, Batch Code
	Caution

## Appendix 1. Electromagnetic Compatibility

Like other electrical medical equipment the system requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC) the system must be installed and operated according to the EMC information provided in this manual.

The system has been designed and tested to comply with IEC 60601-2 requirements for EMC with other devices.

<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>		
The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.		
<b>Emission Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment Guidance</b>
RF emission CISPR 11	Group 1	The system use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emission IEC/EN 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

### Guidance and manufacturer's declaration - electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 test level	Compliance Level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2, 4, 8, 15 kV air	+/- 8 kV contact +/- 2, 4, 8, 15 kV air	If floors are covered with synthetic material the relative humidity shall be least 30%.
Electrical fast transient / burst IEC 61000-4-4	+/- 2 kV for mains supply lines +/- 1 kV for input / output lines	+/- 2 kV mains supply lines N/A	Mains power quality shall be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 0.5, 1 kV line(s) to line(s) +/- 0.5, 1, 2 kV line(s) to earth	+/- 0.5, 1 kV line(s) to line(s), mains supply +/- 0.5, 1, 2 kV line(s) to earth, mains supply	Mains power quality shall be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% $U_t$ (100% dip in $U_t$ ) for 0.5 cycle at 8 angles  0% $U_t$ (100% dip in $U_t$ ) for 1 cycle  70% $U_t$ (30% dip in $U_t$ ) for 25 cycles  0% $U_t$ (100% dip in $U_t$ ) for 5 sec.	100% reduction 0.5 period at 8 angles  100% reduction for 1 period  30% reduction for 25 periods  100% reduction for 5 sec.	Mains power quality shall be that of a typical commercial or hospital environment.  If the use of the system requires continued operation during power mains interruptions the system can be powered by the built in rechargeable battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE:  $U_t$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted Radio Frequency IEC 61000-4-6	3 V RMS 0,15 MHz - 80 MHz  6 V RMS in ISM bands  80% AM at 1 kHz	3 V RMS 0,15 MHz - 80 MHz  6 V RMS in ISM bands  80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the system, including its cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated Radio Frequency IEC 61000-4-3	3 V/m 80-2700 MHz 80% AM at 1 kHz  Additional proximity field test frequencies and levels according to 60601-1-2 table 9	3 V/m 80-2700 MHz 80% AM at 1 kHz  Additional proximity field test frequencies and levels according to 60601-1-2 table 9	Recommended separation distance $d = 1.17 \cdot \sqrt{P}$ $d = 1.17 \cdot \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \cdot \sqrt{P}$ 800 MHz to 2.7 GHz  Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range b).  Interference may occur in the vicinity of equipment marked with the following symbol.  

NOTE 1: At 80 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with a accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey shall be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system shall be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths shall be less than 3 V/m.

### Recommended Separation Distances Between Portable and Mobile RF Communication Equipment and system

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters and the system as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.33\sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.37 m	0.37 m	0.74 m
1	1.17 m	1.17 m	2.33 m
10	3.70 m	3.70 m	7.37 m
100	11.7 m	11.7 m	23.3 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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