
VIVASIGHT

By ETVIEW

Instructions for Use

VivaSight-EB

For use by trained personnel only

EN	Instructions for use - For use by trained personnel only	4
BG	Упътване за употреба - За употреба само от обучен персонал.....	8
CS	Návod k použití - Určeno pouze k použití zaškolenými osobami	12
DA	Brugsvejledning - Må kun anvendes af uddannet personale.....	16
DE	Bedienungsanleitung - Nur zur Verwendung durch geschultes Personal.....	20
EL	Οδηγίες χρήσης - Για χρήση μόνο από εκπαιδευμένο προσωπικό	24
ES	Instrucciones de Uso - para uso exclusivo de personal capacitado.....	28
ET	Kasutusjuhised - Ainult väljaõppinud personalile kasutamiseks	32
FI	Käyttöohjeet - Vain koulutetun henkilöstön käyttöön	36
FR	Mode d'Emploi - Conçu pour du personnel formé uniquement	40
HR	Upute za uporabu - Uporaba je namijenjena samo obučenom osoblju.....	44
HU	Használati utasítás - Kizárólag szakképzett személyzet általi használatra	48
IT	Istruzioni per l'uso - Per l'uso esclusivo da parte di personale addestrato	52
JA	取扱説明書 - 必ず訓練を受けた担当者のみが使用すること	56
KO	사용 지침 - 훈련을 받은 전문인만 사용할 수 있음	60
LT	Naudojimo instrukcijos - Skirtas naudoti tik apmokytiems darbuotojams.....	64
LV	Lietošanas instrukcijas - Driest lietot tikai apmācīts personāls.....	68
NL	Gebruiksaanwijzing - Alleen worden gebruikt door getrainde medewerkers.....	72
NO	Bruksanvisning - Skal kun brukes av opplært personell	76
PL	Instrukcja obsługi - wyłącznie do użytku przez przeszkolony personel.....	80
PT	Instruções de utilização - Para utilização exclusiva de pessoal qualificado.....	84
RO	Instrucțiuni de utilizare - a se utiliza numai de către personal calificat.....	88
RU	Инструкция по эксплуатации - только для обученного персонала.....	92
SK	Návod na použitie - určené na použitie výlučne vyškoleným personálom	96
SL	Navodila za uporabo - le za usposobljeno osebje.....	100
SV	Användarinstruktioner - Uteslutande för användning av utbildad personal.....	104
TR	Kullanım talimatları - sadece eğitimli personelin kullanımına uygundur	108
ZH	使用说明 - 只能由经过培训的人员操作	112

VIVASIGHT-EB

STERILE (ETO) unless package has been opened or damaged.


FOR SINGLE USE ONLY.

DO NOT REUSE.

Federal (US) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

VivaSight-EB is a single-use endobronchial blocker designated for one-lung ventilation. VivaSight-EB is used in combination with VivaSight-SL or equivalent Endotracheal tube. VivaSight-EB features a blue silicone cuff. VivaSight-EB can be deflated and repositioned under visual control provided by VivaSight-SL. The angled distal tip of VivaSight-EB can be adjusted to facilitate navigation to the bronchi.

INDICATIONS FOR USE

VivaSight-EB is used to isolate the left or right lung of a patient in surgery, that requires one lung ventilation or one lung anesthesia.

CONTRAINDICATIONS

VivaSight-EB is contraindicated for use in procedures involving the use of laser equipment in the immediate area of the device.

VivaSight-EB is not recommended for use on patients with serious throat edema or inflammation, hemorrhage, neck vertebra trauma or wet lung.

WARNINGS

- VivaSight-EB is a single use product. Do not reuse: After use, handle and dispose according to local regulations for disposal of infected medical devices. Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.
- Expert clinical judgment should be used in selecting the suitable tracheal tube size for each patient. Use the same size VivaSight-SL as the standard tube size required in combination with VivaSight-EB.
- Make sure the cuff is completely deflated prior to inserting or removing VivaSight-EB.
- Do not attempt to resterilize any part of VivaSight-EB.
- Prior to use, check that the pouch seal is intact and there are no impurities or damage on the product such as rough surfaces, sharp edges or protrusions which may harm the patient.
- Do not use VivaSight-EB if the cuff is damaged. Care must be taken to avoid damaging the cuff during use.

- The cuff pressure should not exceed 20 cm H₂O. Diffusion of a nitrous oxide mixture, oxygen, or air may either increase or decrease the cuff pressure.
- Do not overinflate the cuff. Overinflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or cuff distortion, which may lead to airway blockage.
- Make sure to insert the blocker's tip within the VivaSight-SL or equivalent before connecting the multi joint connector on the blocker to VivaSight-SL or equivalent 15mm connector to avoid pulling the tip back in the multi-joint connector.

CAUTIONS

- VivaSight-EB is intended for use by trained personnel only.
- Use a water-soluble lubricant only. Other lubricants might affect the cuff.
- The use of lidocaine topical aerosols has been associated with the formation of pinholes in PVC cuffs*. To prevent cuff leaks, expert clinical judgment is required when using lidocaine topical aerosols.
* Jayasuriya KD, Watson WF: "P.V.C. cuffs and Lidocaine-based aerosol"; Br J Anaesth. 1981 Dec; 53 (12): 1368.
- Use a cuff pressure gauge to help monitor and adjust the cuff pressure. Inflation of the cuff by "feel" alone or by using a measured volume of air is not recommended as resistance is an unreliable guide during inflation.
- Syringes, stopcocks, or other devices should not be left in the inflation system for an extended period.

STORAGE AND TRANSPORTATION

- Store and transport VivaSight-EB at temperatures between 0° C and 42° C, relative humidity between 10-100% and atmospheric pressure between 80-109 kPa.
- Store in a dry, cool and dark place.

ADVERSE EVENTS

Reported adverse reactions associated with the use of endobronchial blockers are many and diverse. Adverse events associated with the use of VivaSight-EB are the same as those for standard endobronchial blockers. Potential adverse events related to endobronchial blockers include perforating a bronchus or lung parenchyma, causing pneumothorax; mistakenly inflating the bronchial cuff near the tracheal lumen, producing respiratory arrest; developing severe hypoxemia and potential pulmonary edema following continuous suctioning of the nondependent lung by the blocker.

SPECIFICATIONS

- Size: 9 Fr
- Tube length: up to 700 mm
- Compatible for use with VivaSight-SL or equivalent endotracheal tube with inner diameter (mm): 7.0, 7.5, 8.0
- Operating ambient temperature: 10-37 °C (50-98 °F)
- Operating relative humidity: 30-75%
- Operating atmospheric pressure: 80-109 kPa
- Operating altitude: ≤2000 m

BEFORE PLACEMENT

Note: If you are intubating with VivaSight-SL, follow VivaSight-SL instructions for use (IFU).

Do not use if VivaSight-EB's sterile (ETO) package has been damaged or previously opened.

1. Test the cuff for integrity by completely inflating and deflating prior to use.
 - a. For proper inflation, cuff pressure should not exceed 20 cm H₂O.
 - b. For proper deflation, use a 20 cc syringe to completely deflate the cuff.
 - c. Pull the syringe plunger back towards the 20 cc mark.
 - d. Disconnect the syringe from the inflation port without releasing the plunger.

2. Prepare the patient.
3. Lubricate the entire surface of VivaSight-EB's cuff and the interior surface of VivaSight-SL or equivalent.

PLACEMENT

1. Intubate the patient with VivaSight-SL or equivalent.
2. Deflate VivaSight-EB's cuff completely using a 20 cc syringe.
3. Lubricate VivaSight-EB's cuff
4. Keep the VivaSight-EB cuff close, but not in, the multi-joint connector. Connect the VivaSight-EB tube connector to VivaSight-SL or equivalent.
5. Connect the ventilator to the anesthesia circuit on the multi-joint connector.
6. Push VivaSight-EB into VivaSight-SL or equivalent. If resistance is felt during insertion, pull VivaSight-EB backward (avoiding pulling the tip into the multi-joint connector), and then push forward while gently twisting VivaSight-EB.
7. Guide VivaSight-EB's cuff to the target bronchus. The cuff's direction can be changed by twisting the blocker tube left or right.
8. After confirming the position of VivaSight-EB cuff, using VivaSight-SL camera, inflate the cuff to block the target bronchus and lock VivaSight-EB's tube with the twist lock.
9. If required, open the evacuation port following placement to deflate the lung or administer oxygen.

REMOVAL

Deflate the VivaSight-EB cuff completely before removing VivaSight-EB. Follow proper protocol prior to removal.

1. Use a 20 cc syringe to completely deflate the cuff.
2. Pull the syringe plunger back towards the 20 cc mark.
3. Disconnect the syringe from the inflation port without releasing the plunger.
4. Disconnect the VivaSight-EB from the VivaSight-SL or equivalent, and pull it gently backward to extract it from the SL tube.

DISPOSAL

Discard VivaSight-EB in an appropriate biohazard receptacle in accordance with local regulations.

MANUFACTURED FOR:

ETView Ltd.
 Catom 2 street
 Misgav Business Park
 M.P. Misgav 2017900, Israel
www.etview.com

LB00079 REV09

Manufacture Date

Expiry Date

Consult instructions
for use

Do not use if package
is damaged

Do not reuse

Sterilized Using
Ethylene Oxide

US Federal law restricts
this device to sale by
or on the order of a
physician

This product is
not made with phthalates

Temperature for
transportation and
storage

Authorized
Representative in the
European Community

Manufacturer

This product is not made
with natural rubber latex

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